

# Environmental Monitoring by Universal PCR

---

Effective environmental monitoring using a simple assay for measuring and delineating DNA signatures in samples.

Kathy C. Terlesky, Ph.D.  
General Dynamics-Advanced Information Systems  
Charlottesville, VA

# Environmental Monitoring for Target Agents

## NIAID Category A & B Priority Pathogens

### Category A

*Bacillus anthracis*  
*Clostridium botulinum*  
*Yersinia pestis*  
*Francisella tularensis*  
>16 viruses

### Category B

*Burkholderia pseudomallei*  
*Coxiella burnetti*  
*Brucella* species  
*Burkholderia mallei*  
*Rickettsia prowazekii*  
Diarrheagenic *E.coli*  
Pathogenic *Vibrios*  
*Shigella* species  
*Salmonella*  
*Listeria monocytogenes*  
*Campylobacter jejuni*  
*Yersinia enterocolitica*  
*Cryptosporidium parvum*  
*Cyclospora cayatanensis*  
*Giardia lamblia*  
*Entamoeba histolytica*  
*Toxoplasma*  
*Microsporidia*  
+ 11 viruses + 3 toxins


One method of environmental monitoring for homeland defense applications is to monitor the levels of specific threat agents in a target environment.

This requires:

1. The development and validation of multiple agent-specific assays.
2. The application of multiple agent specific assays in each type of environment.

# Environmental Monitoring for Biological Signatures

1. Examine the targeted environment for *multiple environmental signatures* including physical, chemical, biological, and climactic.
2. Analyze data and temporal changes.
3. Develop reliable algorithms and signature processing methods to detect a signal anomaly in the environment.




# **Environmental Monitoring Assay for Biological Signatures**

---

## **Requirements of Assay for Environmental Monitoring of Biological Signatures**

- 1. Simultaneous analysis of multiple (>10) biological signatures in a sample.**  
**One assay can be used in multiple types of environments of known and unknown composition.**
- 2. Detection of target and background populations**
- 3. Detection of genetically-altered agents undetectable by agent-specific assays**
- 4. Cost-effective**

# Universal PCR Methodology




Target gene selected is 16S/18S rRNA.  
(Present in all organisms excluding viruses)

## 1. Universal Amplification

Design primers to a conserved region in target gene.

## 2. Taxon-specific Detection

Design probes to variable regions in target gene.



Region of gene with high sequence conservation.

Region of gene with low sequence conservation.

TAQMAN probes have been developed to detect amplification of target DNA from known classes of organisms.

Bacilli Plant Fungi Proteobacteria

# Universal Amplification

## Primer Validation

Representative Organisms demonstrating a positive response in a SYBR-Green  
Universal PCR Reaction

### *Bacteria-*

*Firmicutes*

*Bacteroidetes*

*Proteobacteria*

### *Archaea*

*Euryarchaeotoa*

### *Eucarya/Fungi*

*Ascomycota*

### *Eucarya/Plant*

*Embryophyta*

### *Eucarya/Animal*

*Chordata*

# Taxon-specific Detection

## Probe Validation

*TAQMAN probes were tested in the PCR assay using DNA with varying target similarity.*

*Data not Available*

-  **Kingdom-level probe**
-  **Family-level probe**
-  **Class-level probe**

# Universal PCR Probe Validation

Annex A from Biological Weapons Convention List (2001)  
and NIAID Category A and B pathogens

*Data not Available*

Only 3 probes (=PCR reactions) required to target the taxonomic class of ALL the bacteria.

In addition, this method will also detect genetically-modified variants of each of these.

Reduces time, cost and complexity of analysis.



# Environmental Monitoring Data

---

Environmental monitoring  
data from air samples in  
Virginia.



*Data not Available*

# Environmental Monitoring Data

---


Amplicons from SYBR Green PCR reaction were cloned and sequenced from 5 separate days.



*Data not Available*


# Environmental Monitoring Instrumentation Developments

Representative melt curves from the air-environmental study.



Broad melt curve from an environmental sample suggesting the presence of multiple amplicon signatures.

New instrumentation available that is specifically suited to analyze complex melt curves.



# An Environmental Monitoring Assay for Multiple Types of Environments

The molecular fingerprinting approach described here can be applied to any sample containing cellular or DNA material.

However, the value of the retrieved data increases in value as more molecular fingerprints are obtained for a sample or environment.

+ signal  +++ signal

## AIR

Background

Fungi

Plant

Targets

*Bacillus*

*Yersinia*

## WATER

Background

Algae

Plant

Targets

*Cryptosporidia*

*Bacillus*

## SOIL

Background

Fungi

Actinobacteria

Targets

*Bacillus*

*Yersinia*

# Environmental Monitoring PCR vs. Agent Monitoring PCR

---

## Environmental Monitoring PCR

One PCR assay can be used in multiple scenarios.

Indicates the presence of DNA from potential warfare agents.

Indicates levels of background DNA.

Can detect DNA from genetically-modified organisms.

Valuable as a control in diagnostic PCR.

PCR Results are evaluated within the context of the sample.

Provides broad taxonomic data on the content of a sample.

Method has not been developed for use for virus detection.

## Agent Monitoring PCR

Requires multiple PCR assays to target multiple threat agents.

Requires continual assay development to target emerging threat agents.

Provides data on presence or absence of specific threat agent.

Can detect and identify viral, bacterial and protozoan targets.

# Summary

---

**You get what you look for.**

**Obtain as many characteristic signatures as possible from a sample.**

**The accumulation of multiple data sources in environmental monitoring can provide a reliable warning of a suspect biological signature.**

**Biological data**

**Meteorological data**

**Chemical data**

**Physical data**