




U.S. GENOMICS



## Pathogen Identification By DNA Barcoding

Mark Nadel



## Outline

- » DirectLinear™ Analysis
- » DARPA *in silico* Test Scenario
- » Technical Framework
- » Results




## Outline

- » **DirectLinear™ Analysis**
- » DARPA *in silico* Test Scenario
- » Technical Framework
- » Results




## DirectLinear™ Analysis




M205XG1 WIAF-269 WIAF-652 WIAF-270

# GeneEngine™ Platform


Sample prep  
and tagging




Microfluidics



Optics



Data capture  
and analysis



GeneEngine™  
Platform



# GeneEngine™ Trace Files





M205XG1 — WIAF-269 WIAF-652  
WIAF-270

## For Further Information

Eugene Y. Chan, Nuno M. Goncalves, Rebecca A. Haeusler, Amie J. Hatch, Jonathan W. Larson, Anthony M. Maletta, Gregory R. Yantz, Eugene D. Carstea, Martin Fuchs, Gordon G. Wong, Steven R. Gullans, Rudolf Gilmanshin,

### **DNA Mapping Using Microfluidic Stretching and Single-Molecule Detection of Fluorescent Site-Specific Tags**

Genome Research (2004), in press.



## Outline

- » DirectLinear™ Analysis
- » **DARPA *in silico* Test Scenario**
- » Technical Framework
- » Results



## ***In silico* Test Scenario**

- » Targets → *Y. pestis*, *V. cholerae*, *B. subtilis*
- » Background → 80 bacteria and yeast species



## Outline

- » DirectLinear™ Analysis
- » DARPA *in silico* Test Scenario
- » **Technical Framework**
- » Results



M205XG1 — WIAF-269 WIAF-652  
WIAF-270


## Underlying Method

- » Fragment DNA to isolate by size
- » Tag with a “motif”
- » Each fragment has a set of possible “perfect” barcodes
  - Determined by thresholding trace files
  - Motif Map + “parameter settings” + starting point → perfect barcode



M205XG1 — WIAP-269 — WIAP-652  
WIAP-270

# GeneEngine™ Simulator





## Basic Test Strategy

- » Check for perfect barcodes of certain fragments of targets
- » Sound the alarm if we see some predetermined number
- » The challenge is to find an enzyme and tag so that each target will have some “identifier” fragment(s) such that:
  - If the organism is present there is a high probability of seeing one of the identifier’s perfect barcodes.
  - If the organism is not present, there is a very low probability of seeing one of the identifier’s perfect barcodes



## Error Model

- » Each fragment has a set of (imperfect) barcodes
  - Perfect sites
  - Single end mismatch (semm) sites.
  - Single interior mismatch (simm) sites.
- » Non-specific binding.

Example: Tag site **ATGC**

Perfect	<b>ATGC</b>
SEMM	<b>TTGC</b>
SIMM	<b>AAGC</b>
Nonspecific	<b>CGAA</b>



## Adding Probabilities

- » Each fragment has a set of (imperfect) barcodes
  - Perfect sites
    - .95 chance of binding at any specific perfect site.
  - Single end mismatch (semm) sites.
    - .04 chance of binding at any specific semm site.
  - Single interior mismatch (simm) sites.
    - .004 chance of binding at any specific simm site.
- » Non-specific binding.
  - .000001 chance of binding at an arbitrary site.



## Key Notions:

- » Tag 2 fragments  $f$  and  $g$
- »  $\text{resemble}(g, f)$  = probability that if interrogation of  $g$  results in a barcode, it will be a perfect barcode of  $f$ 
  - Best toward 0
- »  $\text{robustness}(f) = \text{resemble}(f, f)$ ; probability that if interrogation of  $f$  results in a barcode, it will be a perfect barcode of  $f$ 
  - Best toward 1




## Outline

- » DirectLinear™ Analysis
- » DARPA *in silico* Test Scenario
- » Technical Framework
- » **Results**



## Results of Fragmentation by Asc I





M205XG1 — WIAF-269 — WIAF-652 — WIAF-270

## Identifier Fragment Maps

Fragment at recognition sequence  $gg^{\wedge}cgcg\_cc$   
Tag at agagagag

*Y. pestis* —————++ 214kb

*B. subtilis* —————|————— 295kb

*V. cholerae* ———+——+——+—— 210kb




## Relevant Robustness and Resemble Values

- » Robustness of target identifiers
  - `robustness(pestis_213,776)=0.557.`
  - `robustness(subtilis_295,197)=0.319.`
  - `robustness(cholorae_209,671)=0.307.`
- » Resemble for potential interlopers
  - `resemble(background_218,266, pestis_213,776)= 0`
  - `resemble(cholorae_209,671, pestis_213,776)= 9.68x10-7`
  - `resemble(background_291,301,subtilis_295,197)=2.18x10-9`
  - `resemble(pestis_213,776, cholorae_209,671)=6.42x10-10.`
- » Use efficiency factor of 0.3
  - Only 30% of fragments from original sample result in a barcode.




M205XG1  
WIAF-269  
WIAF-652  
WIAF-270

## ROC Curve for Target *Y. pestis*






## Alternate ROC Curve for Target






M205XG1  
WIAF-269  
WIAF-652  
WIAF-270

## ROC Curve for Target *V. cholerae*





## ROC Curve for Target *B. subtilis*





M205XG1 — WIAF-269 WIAF-652  
WIAF-270

## Acknowledgements

- » DARPA
- » Dr. Bin Feng
- » Dr. Rudolf Gilmanshin
- » John Harris
- » Dr. Hai-Ping Ko
- » Dr. Michael Shia